

兔抗 THOC3 多克隆抗体

中文名称: 兔抗 THOC3 多克隆抗体

英文名称: Anti-THOC3 rabbit polyclonal antibody

别 名: THO complex 3; THO3; hTREX45

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗原: THOC3

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

This gene encodes a component of the nuclear THO transcription elongation complex, which is part of the larger transcription export (TREX) complex that coupl es messenger RNA processing and export. In humans , the transcription export complex is recruited to the 5'-end of messenger RNAs in a splicing- and cap-de pendent manner. Studies of a related complex in mo use suggest that the metazoan transcription export c omplex is involved in cell differentiation and develop ment. A pseudogene of this gene has been defined on chromosome 5.



Applications:	ELISA, WB, IHC
Name of antibody:	THOC3
Immunogen:	Fusion protein of human THOC3
Full name:	THO complex 3
Synonyms:	THO3; hTREX45
SwissProt:	Q96J01
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	39 kDa
WB Positive control:	Jurkat and PC3 cell lysates
WB Recommended dilution:	1000-5000





