

兔抗 THPO 多克隆抗体

中文名称: 兔抗 THPO 多克隆抗体

英文名称: Anti-THPO rabbit polyclonal antibody

别 名: ML; TPO; MGDF; MKCSF; MPLLG; THCYT1

相关类别: 一抗

储 存: 冷冻(-20℃) 避光

宿 主: Rabbit

抗 原: THPO

反应种属: Human, Mouse, Rat

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格



Name of antibody:	THPO
Immunogen:	Fusion protein of human THPO
Full name:	Thrombopoietin
Synonyms:	ML; TPO; MGDF; MKCSF; MPLLG; THCYT1
SwissProt:	P40225
WB Predicted band size:	38 kDa
WB Positive control:	293 cells
WB Recommended dilution:	500-2000

