

## 小鼠抗 POLR2J2 单克隆抗体

中文名称： 小鼠抗 POLR2J2 单克隆抗体

英文名称： Anti-POLR2J2 mouse monoclonal antibody

别 名： HRPB11B; RPB11b1

相关类别： 一抗

储 存： 冷冻（-20℃） 避光

宿 主： Mouse

抗 原： POLR2J2

反应种属： Human, Dog, Rat, Monkey

标 记 物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

#### Background:

This gene is a member of the RNA polymerase II subunit 11 gene family, which includes three genes in a cluster on chromosome 7q22.1 and a pseudogene on chromosome 7p13. The founding member of this family, DNA directed RNA polymerase II polypeptide J, has been shown to encode a subunit of RNA polymerase I, the polymerase responsible for synthesizing messenger RNA in eukaryotes. This locus produces multiple, alternatively spliced transcripts that potentially express isoforms with distinct C-termini compared to DNA directed RNA polymerase II polypeptide J. Most or all variants are spliced to include additional non-coding exons at the 3' end which makes them candidates for nonsense-mediated

	ed decay (NMD). Consequently, it is not known if this locus expresses a protein or proteins in vivo. [provided by RefSeq].
<b>Applications:</b>	WB
<b>Name of antibody:</b>	POLR2J2
<b>Immunogen:</b>	Fusion protein of human POLR2J2
<b>Full name:</b>	polymerase (RNA) II (DNA directed) polypeptide J2 (POLR2J2)
<b>Synonyms:</b>	HRPB11B; RPB11b1
<b>SwissProt:</b>	Q9GZM3
<b>WB Predicted band size:</b>	13 kDa
<b>WB Positive control:</b>	HepG2, Hela, HT-29, A549, COS7, Jurkat, MDCK, PC12, MCF-7 cell lysates
<b>WB Recommended dilution:</b>	1000-5000