

兔抗 TPD52L3 多克隆抗体

- 中文名称: 兔抗 TPD52L3 多克隆抗体
- 英文名称: Anti-TPD52L3 rabbit polyclonal antibody
- 别名: tumor protein D52 like 3; D55; hD55; NYDSP25
- 相关类别: 一抗
- 储存: 冷冻(-20℃)
- 宿 主: Rabbit
- 抗 原: TPD52L3
- 反应种属: Human
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the tumor protein D52-like family of proteins. These proteins are char acterized by an N-terminal coiled-coil motif that is used to form homo- and heteromeric complexes wi th other tumor protein D52-like proteins. The enco ded protein may play a role in spermatogenesis. Al ternative splicing results in multiple transcript varian ts.
Applications:	ELISA, IHC
Name of antibody:	TPD52L3
Immunogen:	Fusion protein of human TPD52L3



Full name:	tumor protein D52 like 3
Synonyms:	D55; hD55; NYDSP25
SwissProt:	Q96J77
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	100-200

