

兔抗 TPH2(Phospho-Ser19) 多克隆抗体

- 中文名称: 兔抗 TPH2(Phospho-Ser19) 多克隆抗体
- 英文名称: Anti-TPH2(Phospho-Ser19) rabbit polyclonal antibody
- 别名: NTPH; ADHD7
- 相关类别: 一抗
- 储存: 冷冻(-20℃)避光
- 宿 主: Rabbit
- 抗原: TPH2(Phospho-Ser19)
- 反应种属: Human Mouse Rat
- 标记物: Unconjugate
- 克隆类型: Unconjugate

技术规格

r t S Background: r e n b f	This gene encodes a member of the pterin-dependent a romatic acid hydroxylase family. The encoded protein ca talyzes the first and rate limiting step in the biosynthesi s of serotonin, an important hormone and neurotransmi tter. The human genome contains two related tryptopha n hydroxylases, one on chromosome 11p15-p14 and on e on chromosome 12q21. This gene is expressed predo minantly in the brain stem. Mutations in this gene may be associated with psychiatric diseases such as bipolar a ffective disorder and major depression. [provided by Ref Seq, Jul 2008]
---	--



Applications:	WB
Name of antibody:	TPH2(Phospho-Ser19)
Immunogen:	Peptide sequence around phosphorylation site of Serine 19(G-F-S(p)-L-D) derived from Human TPH2.
Full name:	tryptophan hydroxylase 2
Synonyms :	NTPH; ADHD7
SwissProt:	Q8IWU9
WB Predicted band size:	56 kDa
WB Positive control:	HepG2 cells lysates
WB Recommended dilution:	500-1000

