

## 小鼠抗 SMARCC1 单克隆抗体

中文名称： 小鼠抗 SMARCC1 单克隆抗体

英文名称： Anti-SMARCC1 mouse monoclonal antibody

别名： SWI/SNF related, matrix associated, actin dependent regulator of chromatin subfamily c member 1; Rsc8; SRG3; SWI3; BAF155; CRACC1

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： SMARCC1

反应种属： Human

标记物： Unconjugate

克隆类型： Mouse Monoclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the SWI/SNF family of proteins, whose members display helicase and ATPase activities and which are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent chromatin remodeling complex SWI/SNF and contains a predicted leucine zipper motif typical of many transcription factors.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	SMARCC1

<b>Immunogen:</b>	Fusion protein of human SMARCC1
<b>Full name:</b>	SWI/SNF related, matrix associated, actin dependent regulator of chromatin subfamily c member 1
<b>Synonyms:</b>	Rsc8; SRG3; SWI3; BAF155; CRACC1
<b>SwissProt:</b>	Q92922
<b>WB Predicted band size:</b>	123 kDa
<b>WB Positive control:</b>	K562, Jurkat and Hela cell lysates
<b>WB Recommended dilution:</b>	1000-5000

