

## 小鼠抗 TBL1X 单克隆抗体

- 中文名称: 小鼠抗 TBL1X 单克隆抗体
- 英文名称: Anti-TBL1X mouse monoclonal antibody
- 别 名: transducin beta like 1X-linked; EBI; TBL1; SMAP55
- 相关类别: 一抗
- 储 存: 冷冻 (-20℃)
- 宿 主: Mouse
- 抗 原: TBL1X
- 反应种属: Human
- 标记物: Unconjugate
- 克隆类型: Mouse Monoclonal

## 技术规格

	The protein encoded by this gene has sequence similarity with members of the WD40 repeat-containing protein family. The WD 40 group is a large family of proteins, which appear to have a r egulatory function. It is believed that the WD40 repeats mediate protein-protein interactions and members of the family are involv
Background:	ed in signal transduction, RNA processing, gene regulation, vesic ular trafficking, cytoskeletal assembly and may play a role in the control of cytotypic differentiation. This encoded protein is found as a subunit in corepressor SMRT (silencing mediator for retinoid and thyroid receptors) complex along with histone deacetylase 3 protein. This gene is located adjacent to the ocular albinism gen



	e and it is thought to be involved in the pathogenesis of the oc ular albinism with late-onset sensorineural deafness phenotype. F our transcript variants encoding two different isoforms have been found for this gene. This gene is highly similar to the Y chromos ome TBL1Y gene.
Applications:	WB, IF
Name of antibody:	TBL1X
Immunogen:	Fusion protein of human TBL1X
Full name:	transducin beta like 1X-linked
Synonyms:	EBI; TBL1; SMAP55
SwissProt:	O60907
WB Predicted band siz e:	62 kDa
WB Positive control:	SW480 cell lysate
WB Recommended dil ution:	1000-5000
IF positive control:	Hela cells
IF Recommend dilution	50-100





