

## 小鼠抗 ZNF670 单克隆抗体

中文名称：小鼠抗 ZNF670 单克隆抗体

英文名称：Anti-ZNF670 mouse monoclonal antibody

别名：zinc finger protein 670

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：ZNF670

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

### 技术规格

<b>Background:</b>	May be involved in transcriptional regulation.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	ZNF670
<b>Immunogen:</b>	Fusion protein of human ZNF670
<b>Full name:</b>	zinc finger protein 670
<b>SwissProt:</b>	Q9BS34
<b>WB Predicted band size:</b>	45 KD
<b>WB Positive control:</b>	HepG2 cell lysate
<b>WB Recommended dilution:</b>	500-2000