

小鼠抗 ZSCAN4 单克隆抗体

中文名称: 小鼠抗 ZSCAN4 单克隆抗体

英文名称: Anti-ZSCAN4 mouse monoclonal antibody

别名: ZNF494

抗原: ZSCAN4

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human , Dog

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	The ZSCAN4 gene encodes a protein involved in telomere maintenance and with a key role in the critical feature of mouse embryonic stem (ES) cells, namely, defying cellular senescence and maintaining normal karyotype for many cell divisions in culture (Zalzman et al., 2010 [PubMed 20336070]).[supplied by OMIM, May
Applications:	WB
Name of antibody:	ZSCAN4
Immunogen:	Fusion protein of human ZSCAN4
Full name:	zinc finger and SCAN domain containing 4 (ZSCAN4)
Synonyms:	ZNF494
SwissProt:	Q8NAM6
WB Predicted band size:	49 kDa



WB Positive control:	MDCK cell lysate
WB Recommended dilution:	1000-5000