

## SHC3 抗原（重组蛋白）

中文名称： SHC3 抗原（重组蛋白）

英文名称： SHC3 Antigen (Recombinant Protein)

别名： SHC adaptor protein 3; RAI; NSHC; SHCC; N-Shc

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 401-594 amino acids of human SHC3

### 技术规格

<b>Full name:</b>	SHC adaptor protein 3
<b>Synonyms:</b>	RAI; NSHC; SHCC; N-Shc
<b>Swissprot:</b>	Q92529
<b>Gene Accession:</b>	BC026314
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Signaling adapter that couples activated growth factor receptors to signaling pathway in neurons. Involved in the signal transduction pathways of neurotrophin-activated Trk receptors in cortical neurons.