

SKOV3 细胞裂解液

中文名称: SKOV3 细胞裂解液

英文名称: SKOV3 Cell Lysate

别名: SKOV3 细胞裂解液

相关类别: 抗原

储存: 冷冻 (-20℃) 避光

概述

SKOV3 cell lysates from human ovarian cancer cells were cultured in the media containing 10% FBS.

技术规格

Backg round	Cells were lysed in lysis buffer (50mM Tris(PH7.4), 150mM NaCl, 1% NP-40, 1% SD S). 100 µg of cellular protein was mixed with 4 times volume of 5X SDS-PAGE loading buffer.
------------------------	---