

SLC35E2A 抗原（重组蛋白）

中文名称： SLC35E2A 抗原（重组蛋白）

英文名称： SLC35E2A Antigen (Recombinant Protein)

别名： solute carrier family 35 member E2A; SLC35E2

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 187-266 amino acids of human SLC35E2A

技术规格

Full name:	solute carrier family 35 member E2A
Synonyms:	SLC35E2
Swissprot:	P0CK97
Gene Accession:	BC112289
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Putative transporter.

mlbio 梅联生物
Good elisakit producers