

## SMYD3 抗原（重组蛋白）

中文名称： SMYD3 抗原（重组蛋白）

英文名称： SMYD3 Antigen (Recombinant Protein)

别 名： SET and MYND domain containing 3; KMT3E; ZMYND1; ZNFN3A1; bA74P14.1

储 存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to C terminal 228 amino acids of human SMYD3

### 技术规格

<b>Full name:</b>	SET and MYND domain containing 3
<b>Synonyms:</b>	KMT3E; ZMYND1; ZNFN3A1; bA74P14.1
<b>Swissprot:</b>	Q9H7B4
<b>Gene Accession:</b>	BC031010
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a histone methyltransferase which functions in RNA polymerase II complexes by an interaction with a specific RNA helicase. Multiple transcript variants encoding different isoforms have been found for this gene.