

## TBC1D2 抗原（重组蛋白）

中文名称：TBC1D2 抗原（重组蛋白）

英文名称： TBC1D2 Antigen (Recombinant Protein)

别名： TBC1 domain family member 2; PARIS1; PARIS-1; TBC1D2A

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 729-928 amino acids of human TBC1D2

### 技术规格

<b>Full name:</b>	TBC1 domain family member 2
<b>Synonyms:</b>	PARIS1; PARIS-1; TBC1D2A
<b>Swissprot:</b>	Q9BYX2
<b>Gene Accession:</b>	BC028918
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Acts as GTPase-activating protein for RAB7A. Signal effector acting as a linker between RAC1 and RAB7A, leading to RAB7A inactivation and subsequent inhibition of cadherin degradation and reduced cell-cell adhesion.