

TCF25 抗原(重组蛋白)

中文名称: TCF25 抗原(重组蛋白)

英文名称: TCF25 Antigen (Recombinant Protein)

别 名: transcription factor 25; Hulp1; NULP1; FKSG26; PRO2620; hKIAA1049

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 477-676 amino acids of human TCF25

技术规格

Full name:	transcription factor 25
Synonyms:	Hulp1; NULP1; FKSG26; PRO2620; hKIAA1049
Swissprot:	Q9BQ70
Gene Accession:	BC009349
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May play a role in cell death control. Acts as a transcriptional
	repressor. Has been shown to repress transcription of SRF in v
	itro and so may play a role in heart development.