

TGM4 抗原(重组蛋白)

- 中文名称: TGM4 抗原 (重组蛋白)
- 英文名称: TGM4 Antigen (Recombinant Protein)
- 别名: transglutaminase 4; TGP; hTGP
- 储存: 冷冻(-20℃)
- 相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 485-684 amino acids of human TGM4

技术规格

Full name:	transglutaminase 4
Synonyms:	TGP; hTGP
Swissprot:	P49221
Gene Accession:	BC007003
Purity:	>85%, as determined by Coomassie blue stained SD S-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Associated with the mammalian reproductive process. Catalyzes the cross-linking of proteins and the conju gation of polyamines to specific proteins in the semi nal tract.