

## THAP12 抗原（重组蛋白）

中文名称： THAP12 抗原（重组蛋白）

英文名称： THAP12 Antigen (Recombinant Protein)

别 名： THAP domain containing 12; DAP4; THAP0; PRKRIR; P52rIPK

储 存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 612-761 amino acids of human THAP12

### 技术规格

<b>Full name:</b>	THAP domain containing 12
<b>Synonyms:</b>	DAP4; THAP0; PRKRIR; P52rIPK
<b>Swissprot:</b>	O43422
<b>Gene Accession:</b>	BC021992
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Upstream regulator of interferon-induced serine/threonine protein kinase R (PKR). May block the PKR-inhibitory function of DNAJC3, resulting in restoration of kinase activity and suppression of cell growth.