

XPO6 抗原(重组蛋白)

中文名称: XPO6 抗原(重组蛋白)

英文名称: XPO6 Antigen (Recombinant Protein)

别 名: exportin 6; EXP6; RANBP20

相关类别: 抗原

储 存: 冷冻(-20℃)

概述

Fusion protein corresponding to a region derived from 926-1125 amino acids of human XPO6 $\,$

技术规格

Full name:	exportin 6
Synonyms:	EXP6; RANBP20
Swissprot:	Q96QU8
Gene Accession:	BC014071
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the importin-beta family. Members of this family are regulated by the GTPase Ran to mediate transport of cargo across the nuclear envelope. This protein has been shown to mediate nuclear export of profilin-actin complexe s. A pseudogene of this gene is located on the long arm of chromo



some 14. Alternative splicing results in multiple transcript variants th at encode different protein isoforms.