

YY1AP1 抗原（重组蛋白）

中文名称： YY1AP1 抗原（重组蛋白）

英文名称： YY1AP1 Antigen (Recombinant Protein)

别名： YY1 associated protein 1; YAP; HCCA1; HCCA2; YY1AP

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to C terminal 200 amino acids of human YY1AP1

技术规格

Full name:	YY1 associated protein 1
Synonyms:	YAP; HCCA1; HCCA2; YY1AP
Swissprot:	Q9H869
Gene Accession:	BC009358
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The encoded gene product presumably interacts with YY1 protein; however, its exact function is not known. Alternative splicing results in multiple transcript variants encoding different isoforms.