

ZFP36 抗原（重组蛋白）

中文名称： ZFP36 抗原（重组蛋白）

英文名称： ZFP36 Antigen (Recombinant Protein)

别名： ZFP36 ring finger protein; TTP; GOS24; GOS24; TIS11; NUP475; zfp-36; RNF162A

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 127-326 amino acids of human ZFP36

技术规格

Full name:	ZFP36 ring finger protein
Synonyms:	TTP; GOS24; GOS24; TIS11; NUP475; zfp-36; RNF162A
Swissprot:	P26651
Gene Accession:	BC009693
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	mRNA-binding protein involved in post-transcriptional regulation of A U-rich element (ARE)-containing mRNAs. Acts by specifically binding ARE-containing mRNAs and promoting their degradation. Recruits deadenylase CNOT7 (and probably the CCR4-NOT complex) via association with CNOT1. Plays a key role in the post-transcriptional regulation of tumor necrosis factor (TNF). Plays a key role in the post-transcriptional regulation of tumor necrosis factor (TNF).

mlbio 梅联生物
Good elisakit producers