

鸡 lgY(lgG)(H+L)

- 中文名称: 鸡 lgY(lgG) (H+L)
- 英文名称: Chicken IgY(IgG) (H+L)
- 相关类别: IgG
- 储 存: 冷冻 (-20℃) 避光

技术规格

Background	Immunoglobulin G (IgG) is one of the most abundant proteins i n serum with normal levels between 8-17 mg/mL in adult bloo d. IgG is important for our defense against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions, to bind to the pathogen that elicited the response a nd to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombi nation and the number of specificities in an individual at a give n time point is estimated to be 1011 variants.
Applications	Optimal tested applications and working dilutions must be dete rmined by the end user.