

兔抗 HELT 多克隆抗体

中文名称: 兔抗 HELT 多克隆抗体

英文名称: Anti-HELT rabbit polyclonal antibody

别名: helt bHLH transcription factor; Mgn; HESL; HCM1228; bHLHb44

储存: 冷冻 (-20℃)

抗原: HELT

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Transcriptional repressor which binds preferentially to the canonical E box sequence 5'-CACGCG-3'.
Applications:	ELISA, IHC
Name of antibody:	HELT
Immunogen:	Synthetic peptide of human HELT
Full name:	helt bHLH transcription factor
Synonyms:	Mgn; HESL; HCM1228; bHLHb44
SwissProt:	A6NFD8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human tonsil
IHC Recommend dilution:	30-150



