

兔抗 VRK3 多克隆抗体

中文名称：兔抗 VRK3 多克隆抗体

英文名称：Anti-VRK3 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）避光

宿 主：Rabbit

抗 原：VRK3

反应种属：Human

标 记 物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:

This gene encodes a member of the vaccinia-related kinase (VRK) family of serine/threonine protein kinases. In both human and mouse, this gene has substitutions at several residues within the ATP binding motifs that in other kinases have been shown to be required for catalysis. In vitro assays indicate the protein lacks phosphorylation activity. The protein, however, likely retains its substrate binding capability. This gene is widely expressed in human tissues and its protein localizes to the nucleus. Alternative splicing results in multiple transcripts encoding different isoforms.

Applications:	WB
Name of antibody:	VRK3
Immunogen:	Synthesized peptide derived from Internal of human VRK3.
Full name:	vaccinia related kinase 3
Synonyms :	0
SwissProt:	Q8IV63
WB Predicted band size:	53 kDa
WB Positive control:	COLO205 cells lysate
WB Recommended dilution:	500-3000

