

兔抗 WNT10B 多克隆抗体

中文名称： 兔抗 WNT10B 多克隆抗体

英文名称： Anti-WNT10B rabbit polyclonal antibody

别 名： SHFM6; WNT-12

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： WNT10B

反应种属： Human, Mouse

标 记 物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It may be involved in breast cancer, and its protein signaling is likely a molecular switch that governs adipogenesis. This protein is 96 % identical to the mouse Wnt10b protein at the amino acid level.
Applications:	ELISA, WB, IHC
Name of antibody:	WNT10B
Immunogen:	Synthetic peptide of human WNT10B
Full name:	wingless-type MMTV integration site family, member 10

	B
Synonyms :	SHFM6; WNT-12
SwissProt:	O00744
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human brain
IHC Recommend dilution:	10-50
WB Predicted band size:	43 kDa
WB Positive control:	Hela cell
WB Recommended dilution:	500-2000

