

兔抗 XAB2 多克隆抗体

中文名称: 兔抗 XAB2 多克隆抗体

英文名称: Anti-XAB2 rabbit polyclonal antibody

别名: HCNP; HCRN; SYF1; NTC90

储存: 冷冻 (-20℃)

抗原: XAB2

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Involved in transcription-coupled repair (TCR), transcript ion and pre-mRNA splicing.
Applications:	ELISA, IHC
Name of antibody:	XAB2
Immunogen:	Fusion protein of human XAB2
Full name:	XPA binding protein 2
Synonyms :	HCNP; HCRN; SYF1; NTC90
SwissProt:	Q9HCS7
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human prostate cancer
IHC Recommend dilution:	20-100

