

兔抗 ZDHHC20 多克隆抗体

中文名称: 兔抗 ZDHHC20 多克隆抗体

英文名称: Anti-ZDHHC20 rabbit polyclonal antibody

别名: DHHC-20; 4933421L13Rik

储存: 冷冻 (-20℃)

抗原: ZDHHC20

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

Zinc-finger proteins contain DNA-binding domains and have a wide variety of functions, most of which encompass some form of transcriptional activation or repression. ZDHHC20 (zinc finger, DHHC-type containing 20) is a 365 amino acid multi-pass membrane protein that contains one DHHC-type zinc finger and is thought to function as a palmitoyltransferase, catalyzing the transformation of palmitoyl-CoA and a cysteine-conjugated protein to an S-palmitoyl protein and free CoA. In response to DNA damage, ZDHHC20, which exists as multiple alternatively spliced isoforms, may be phosphorylated by ATM or ATR. The gene encoding ZDHHC20 maps to human c

	hromosome 13, which houses over 400 genes, such as BRCA2 and RB1, and comprises nearly 4% of the human genome.
Applications:	ELISA, IHC
Name of antibody:	ZDHHC20
Immunogen:	Fusion protein of human ZDHHC20
Full name:	zinc finger, DHHC-type containing 20
Synonyms :	DHHC-20; 4933421L13Rik
SwissProt:	Q5W0Z9
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human esophagus cancer
IHC Recommend dilution:	25-100

