

兔抗 ZFP28 多克隆抗体

中文名称: 兔抗 ZFP28 多克隆抗体

英文名称: Anti-ZFP28 rabbit polyclonal antibody

别名: mkr5

储存: 冷冻 (-20℃)

抗原: ZFP28

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	May be involved in transcriptional regulation. May have a role in embryonic development.
Applications:	ELISA, WB
Name of antibody:	ZFP28
Immunogen:	Fusion protein of human ZFP28
Full name:	ZFP28 zinc finger protein
Synonyms :	mkr5
SwissProt:	Q8NHY6
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	99 kDa
WB Positive control:	293T cell
WB Recommended dilution:	200-1000

