

兔抗 ZNF787 多克隆抗体

中文名称：兔抗 ZNF787 多克隆抗体

英文名称：Anti-ZNF787 rabbit polyclonal antibody

级 别：TIP20

抗 原：ZNF787

储 存：冷冻（-20℃）

宿 主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in transcriptional regulation.
Applications:	ELISA, IHC
Name of antibody:	ZNF787
Immunogen:	Synthetic peptide of human ZNF787
Full name:	zinc finger protein 787
Synonyms:	TIP20
SwissProt:	Q6DD87
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human esophagus cancer
IHC Recommend dilution:	30-150

