

兔抗 GUSB 多克隆抗体

中文名称：兔抗 GUSB 多克隆抗体

英文名称：Anti-GUSB rabbit polyclonal antibody

别名：glucuronidase beta; BG; MPS7

抗原：GUSB

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a hydrolase that degrades glycosaminoglycans, including heparan sulfate, dermatan sulfate, and chondroitin-4,6-sulfate. The enzyme forms a homotetramer that is localized to the lysosome. Mutations in this gene result in mucopolysaccharidosis type VII. Alternative splicing results in multiple transcript variants. There are many pseudogenes of this locus in the human genome.
Applications:	ELISA, WB, IHC
Name of antibody:	GUSB
Immunogen:	Fusion protein of human GUSB

Full name:	glucuronidase beta
Synonyms:	BG; MPS7
SwissProt:	P08236
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	40-200
WB Predicted band size:	75 kDa
WB Positive control:	K562 cell lysate
WB Recommended dilution:	500-2000



