

兔抗 GTF2H4 多克隆抗体

- 中文名称：兔抗 GTF2H4 多克隆抗体
- 英文名称：Anti-GTF2H4 rabbit polyclonal antibody
- 别名：general transcription factor IIH subunit 4; P52; TFB2; TFIIH
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：GTF2H4
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	Component of the core-TFIIH basal transcription factor involved in nucleotide excision repair (NER) of DNA and, when complexed to CAK, in RNA transcription by RNA polymerase II.
Applications:	ELISA, IHC
Name of antibody:	GTF2H4
Immunogen:	Synthetic peptide of human GTF2H4
Full name:	general transcription factor IIH subunit 4
Synonyms:	P52; TFB2; TFIIH
SwissProt:	Q92759
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and Human esophagus cancer

IHC Recommend dilution:	40-200
WB Predicted band size:	52 kDa
WB Positive control:	Human heart tissue lysate
WB Recommended dilution:	200-1000

