

兔抗 RBM12 多克隆抗体

中文名称:	兔抗 RBM12 多克隆抗体
英文名称:	Anti-RBM12 rabbit polyclonal antibody
别名:	RNA binding motif protein 12; SWAN; SCZD19; HRIHFB2091
储存:	冷冻 (-20℃)
抗原:	RBM12
宿主:	Rabbit
反应种属:	Human, Mouse
相关类别:	一抗
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	This gene encodes a protein that contains several RNA-binding motifs, potential transmembrane domains, and proline-rich regions. This gene and the gene for copine I overlap at map location 20q11.21. Alternative splicing in the 5' UTR results in four transcript variants. All variants encode the same protein.
Applications:	ELISA, IHC
Name of antibody:	RBM12
Immunogen:	Fusion protein of human RBM12
Full name:	RNA binding motif protein 12
Synonyms:	SWAN; SCZD19; HRIHFB2091

SwissProt:	Q9NTZ6
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and Human prostate cancer
IHC Recommend dilution:	50-300

