

兔抗 RTCB 多克隆抗体

中文名称：兔抗 RTCB 多克隆抗体

英文名称：Anti-RTCB rabbit polyclonal antibody

别名：RNA 2',3'-cyclic phosphate and 5'-OH ligase; FAAP; HSPC117; C22orf28; DJ149A16.6

相关类别：一抗

抗原：RTCB

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Catalytic subunit of the tRNA-splicing ligase complex that acts by directly joining spliced tRNA halves to mature-sized tRNAs by incorporating the precursor-derived splice junction phosphate into the mature tRNA as a canonical 3',5'-phosphodiester. May act as an RNA ligase with broad substrate specificity, and may function toward other RNAs.
Applications:	ELISA, IHC
Name of antibody:	RTCB
Immunogen:	Fusion protein of human RTCB
Full name:	RNA 2',3'-cyclic phosphate and 5'-OH ligase
Synonyms:	FAAP; HSPC117; C22orf28; DJ149A16.6

SwissProt:	Q9Y3I0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and human prostate cancer
IHC Recommend dilution:	40-200



