

兔抗 SNAPC1 多克隆抗体

中文名称：兔抗 SNAPC1 多克隆抗体

英文名称：Anti-SNAPC1 rabbit polyclonal antibody

别名：small nuclear RNA activating complex polypeptide 1; SNAP43; PTFgamma

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SNAPC1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Part of the SNAPc complex required for the transcription of both RNA polymerase II and III small-nuclear RNA genes. Binds to the proximal sequence element (PSE), a non-TATA-box basal promoter element common to these 2 types of genes. Recruits TBP and BRF2 to the U6 snRNA TATA box.
Applications:	ELISA, IHC
Name of antibody:	SNAPC1
Immunogen:	Fusion protein of human SNAPC1
Full name:	small nuclear RNA activating complex polypeptide 1
Synonyms:	SNAP43; PTFgamma
SwissProt:	Q16533

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and Human cervical cancer
IHC Recommend dilution:	50-300



