

兔抗 NDUFB5 多克隆抗体

中文名称：兔抗 NDUFB5 多克隆抗体

英文名称：Anti-NDUFB5 rabbit polyclonal antibody

别名：SGDH; CISGDH

相关类别：一抗

抗原：NDUFB5

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a subunit of the mitochondrial NADH:ubiquinone oxidoreductase (complex I). Mammalian complex I is composed of 45 different subunits. It locates at the mitochondrial inner membrane. This protein has NADH dehydrogenase activity and oxidoreductase activity. It transfers electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone. Three transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	NDUFB5

Immunogen:	Full length fusion protein
Full name:	NADH:ubiquinone oxidoreductase subunit B5
Synonyms:	SGDH; CISGDH
SwissProt:	O43674
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	20-100

