

兔抗 TBL1X 多克隆抗体

中文名称: 兔抗 TBL1X 多克隆抗体

英文名称: Anti-TBL1X rabbit polyclonal antibody

别名: EBI; TBL1; SMAP55

储存: 冷冻 (-20℃)

抗原: TBL1X

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene has sequence similarity with members of the WD40 repeat-containing protein family. The WD40 group is a large family of proteins, which appear to have a regulatory function. It is believed that the WD40 repeats mediate protein-protein interactions and members of the family are involved in signal transduction, RNA processing, gene regulation, vesicular trafficking, cytoskeletal assembly and may play a role in the control of cytotypic differentiation.
Applications:	ELISA, WB
Name of antibody:	TBL1X

Immunogen:	Fusion protein of human TBL1X
Full name:	transducin (beta)-like 1X-linked
Synonyms :	EBI; TBL1; SMAP55
SwissProt:	O60907
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	62 kDa
WB Positive control:	Human lung cancer tissue
WB Recommended dilution:	200-1000

