

兔抗 VILL 多克隆抗体

中文名称：兔抗 VILL 多克隆抗体

英文名称：Anti-VILL rabbit polyclonal antibody

相关类别：一抗

抗原：VILL

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|------------------------------------|---|
| Background: | The protein encoded by this gene belongs to the villin/gelsolin family. It contains 6 gelsolin-like repeats and a headpiece domain. It may play a role in actin-bundling. |
| Applications: | ELISA, IHC |
| Name of antibody: | VILL |
| Immunogen: | Fusion protein of human VILL |
| Full name: | villin-like |
| SwissProt: | O15195 |
| ELISA Recommended dilution: | 5000-10000 |
| IHC positive control: | Human brain and human lung cancer |
| IHC Recommend dilution: | 50-350 |

