

兔抗 POLR1E 多克隆抗体

中文名称：兔抗 POLR1E 多克隆抗体

英文名称：Anti-POLR1E rabbit polyclonal antibody

别名：RNA polymerase I subunit E; PAF53; PRAF1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：POLR1E

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. Component of RNA polymerase I which synthesizes ribosomal RNA precursors. Appears to be involved in the formation of the initiation complex at the promoter by mediating the interaction between Pol I and U BTF/UBF (By similarity).
Applications:	ELISA, IHC
Name of antibody:	POLR1E
Immunogen:	Fusion protein of human POLR1E
Full name:	RNA polymerase I subunit E

Synonyms:	PAF53; PRAF1
SwissProt:	Q9GZS1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human tonsil
IHC Recommend dilution:	40-200

