

兔抗 PPP3R1 多克隆抗体

中文名称：兔抗 PPP3R1 多克隆抗体

英文名称：Anti-PPP3R1 rabbit polyclonal antibody

别名：protein phosphatase 3 regulatory subunit B, alpha; CNB; CNB1; CALNB1

相关类别：一抗

储存：冷冻（-20℃）

抗原：PPP3R1

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Regulatory subunit of calcineurin, a calcium-dependent, calmodulin stimulated protein phosphatase. Confers calcium sensitivity.
Applications:	ELISA, IHC
Name of antibody:	PPP3R1
Immunogen:	Full length fusion protein
Full name:	protein phosphatase 3 regulatory subunit B, alpha
Synonyms:	CNB; CNB1; CALNB1
SwissProt:	P63098
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	50-300

