

兔抗 PRPF3 多克隆抗体

中文名称:	兔抗 PRPF3 多克隆抗体
英文名称:	Anti-PRPF3 rabbit polyclonal antibody
别名:	pre-mRNA processing factor 3; PRP3; RP18; HPRP3; Prp3p; HPRP3P; SNRNP90
抗原:	PRPF3
储存:	冷冻 (-20℃)
宿主:	Rabbit
反应种属:	Human, Mouse
相关类别:	一抗
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	The removal of introns from nuclear pre-mRNAs occurs on complexes called spliceosomes, which are made up of 4 small nuclear ribonucleoprotein (snRNP) particles and an undefined number of transiently associated splicing factors. This gene product is one of several proteins that associate with U4 and U6 snRNPs. Mutations in this gene are associated with retinitis pigmentosa-18.
Applications:	ELISA, WB, IHC
Name of antibody:	PRPF3
Immunogen:	Fusion protein of human PRPF3
Full name:	pre-mRNA processing factor 3

Synonyms:	PRP3; RP18; HPRP3; Prp3p; HPRP3P; SNRNP90
SwissProt:	O43395
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and Human colorectal cancer
IHC Recommend dilution:	20-100
WB Predicted band size:	78 kDa
WB Positive control:	HepG2 cell lysate
WB Recommended dilution:	1000-5000



