

兔抗 PRPF40A 多克隆抗体

中文名称：兔抗 PRPF40A 多克隆抗体

英文名称：Anti-PRPF40A rabbit polyclonal antibody

别名：pre-mRNA processing factor 40 homolog A; HYPA; FBP11; FLAF1; FNBP3; HIP10; Prp40; FBP-11; HIP-10; NY-REN-6

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PRPF40A

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Binds to WASL/N-WASP and suppresses its translocation from the nucleus to the cytoplasm, thereby inhibiting its cytoplasmic function (By similarity). Plays a role in the regulation of cell morphology and cytoskeletal organization. Required in the control of cell shape and migration. May play a role in cytokinesis. May be involved in pre-mRNA splicing.
Applications:	ELISA, IHC
Name of antibody:	PRPF40A
Immunogen:	Fusion protein of human PRPF40A

Full name:	pre-mRNA processing factor 40 homolog A
Synonyms:	HYPA; FBP11; FLAF1; FNBP3; HIP10; Prp40; FBP-11; HIP-10; NY-REN-6
SwissProt:	O75400
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human thyroid cancer
IHC Recommend dilution:	50-200



