

兔抗 PTGDR2 多克隆抗体

中文名称: 兔抗 PTGDR2 多克隆抗体

英文名称: Anti-PTGDR2 rabbit polyclonal antibody

别名: prostaglandin D2 receptor 2; DP2; DL1R; CD294; CRTH2; GPR44

储存: 冷冻 (-20℃)

抗原: PTGDR2

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a G-protein-coupled receptor that is preferentially expressed in CD4+ effector T helper 2 (Th2) cells. This protein is a prostaglandin D2 receptor that mediates the pro-inflammatory chemotaxis of eosinophils, basophils, and Th2 lymphocytes generated during allergic inflammation. Single nucleotide polymorphisms in the 3' UTR of this gene have been associated with asthma susceptibility.
Applications:	ELISA, IHC
Name of antibody:	PTGDR2
Immunogen:	Synthetic peptide of human PTGDR2

Full name:	prostaglandin D2 receptor 2
Synonyms:	DP2; DL1R; CD294; CRTH2; GPR44
SwissProt:	Q9Y5Y4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human ovarian cancer
IHC Recommend dilution:	30-150



