

兔抗 ABLIM1 多克隆抗体

中文名称: 兔抗 ABLIM1 多克隆抗体

英文名称: Anti-ABLIM1 rabbit polyclonal antibody

别 名: ABLIM; LIMAB1; LIMATIN; abLIM-1

抗 原: ABLIM1

储 存: 冷冻(-20℃)

宿 主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a cytoskeletal LIM protein that binds to actin filaments via a domain that is homologous to erythrocyte dematin. LIM domains, found in over 60 proteins, play key roles in the regulation of developmental pathways. LIM domains also function as protein-binding interfaces, mediating specific protein-protein interactions. The protein encoded by this gene could mediate such interactions between actin filaments and cytoplasmic targets. Alternatively spliced transcript variants encoding different isoforms have been identified.
Applications:	ELISA, WB, IHC



Name of antibody:	ABLIM1
Immunogen:	Fusion protein of human ABLIM1
Full name:	Actin binding LIM protein 1
Synonyms:	ABLIM; LIMAB1; LIMATIN; abLIM-1
SwissProt:	O14639
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	15-50
WB Predicted band size:	88 kDa
WB Positive control:	Mouse heart tissue
WB Recommended dilution:	500-2000





