

兔抗 ANO7 多克隆抗体

中文名称：兔抗 ANO7 多克隆抗体

英文名称：Anti-ANO7 rabbit polyclonal antibody

别 名：NGEP; DTMPP; IPCA5; D-TMPP; IPCA-5; PCANAP5; TMEM16G; PCANAP5L

相关类别：一抗

抗 原：ANO7

储 存：冷冻（-20℃）避光

宿 主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This prostate-specific gene encodes a cytoplasmic protein, as well as a polytopic membrane protein which may serve as a target in prostate cancer diagnosis and immunotherapy. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	WB, IHC, IF
Name of antibody:	ANO7
Immunogen:	Synthesized peptide derived from C-terminal of human TM16G.
Full name:	anoctamin 7
Synonyms :	NGEP; DTMPP; IPCA5; D-TMPP; IPCA-5; PCANAP5; TMEM16G; PCANAP5L

SwissProt:	Q6IWH7
IHC positive control:	Human prostate carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	106 kDa
WB Positive control:	LOVO cells lysate
WB Recommended dilution:	500-3000
IF positive control:	HepG2 cells
IF Recommend dilution:	100-500





