

兔抗 BIRC3 多克隆抗体

中文名称：兔抗 BIRC3 多克隆抗体

英文名称：Anti-BIRC3 rabbit polyclonal antibody

别名：baculoviral IAP repeat containing 3; AIP1; API2; MIHC; CIAP2; HAIP1; HIAP1; IAP-1; MALT2; RNF49; c-IAP2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BIRC3

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the IAP family of proteins that inhibit apoptosis by binding to tumor necrosis factor receptor-associated factors TRAF1 and TRAF2, probably by interfering with activation of ICE-like proteases. The encoded protein inhibits apoptosis induced by serum deprivation but does not affect apoptosis resulting from exposure to menadione, a potent inducer of free radicals. It contains 3 baculovirus IAP repeats and a ring finger domain. Transcript variants encoding the same isoform have been identified.
Applications:	ELISA, WB, IHC

Name of antibody:	BIRC3
Immunogen:	Fusion protein of human BIRC3
Full name:	baculoviral IAP repeat containing 3
Synonyms:	AIP1; API2; MIHC; CIAP2; HAIP1; HIAP1; IAP-1; MALT2; RNF 49; c-IAP2
SwissProt:	Q13489
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human tonsil
IHC Recommend dilution:	100-200
WB Predicted band size:	68 kDa
WB Positive control:	Raji cell lysate
WB Recommended dilution:	500-2000





