

兔抗 BRS3 多克隆抗体

中文名称：兔抗 BRS3 多克隆抗体

英文名称：Anti-BRS3 rabbit polyclonal antibody

别名：BB3

相关类别：一抗

抗原：BRS3

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a G protein-coupled membrane receptor that binds bombesin-like peptides. This binding results in activation of a phosphatidylinositol-calcium second messenger system, with physiological effects including regulation of metabolic rate, glucose metabolism, and hypertension.
Applications:	ELISA, WB
Name of antibody:	BRS3
Immunogen:	Synthetic peptide of human BRS3
Full name:	bombesin like receptor 3
Synonyms:	BB3
SwissProt:	P32247

ELISA Recommended dilution:	5000-10000
WB Predicted band size:	44 kDa
WB Positive control:	A549 cell lysate
WB Recommended dilution:	200-1000

