

兔抗 BTG3 多克隆抗体

中文名称：兔抗 BTG3 多克隆抗体

英文名称：Anti-BTG3 rabbit polyclonal antibody

别 名：ANA, TOB5, TOFA, TOB55

储 存：冷冻（-20℃）

抗 原：BTG3

宿 主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the BTG/Tob family. This family has structurally related proteins that appear to have antiproliferative properties. This encoded protein might play a role in neurogenesis in the central nervous system. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	BTG3
Immunogen:	Fusion protein of human BTG3
Full name:	BTG anti-proliferation factor 3
Synonyms:	ANA; TOB5; TOFA; APRO4; TOB55; ANA/BTG3
SwissProt:	Q14201

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	100-200
WB Predicted band size:	29 kDa
WB Positive control:	Mouse testis tissue lysate
WB Recommended dilution:	1000-5000



