

兔抗 CATSPERE 多克隆抗体

中文名称: 兔抗 CATSPERE 多克隆抗体

英文名称: Anti-CATSPERE rabbit polyclonal antibody

别名: catsper channel auxiliary subunit epsilon; C1orf10

储存: 冷冻 (-20℃)

抗原: CATSPERE

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Auxiliary component of the CatSper complex, a complex involved in sperm cell hyperactivation. Sperm cell hyperactivation is needed for sperm motility which is essential late in the preparation of sperm for fertilization.
Applications:	ELISA, IHC
Name of antibody:	CATSPERE
Immunogen:	Fusion protein of human CATSPERE
Full name:	catsper channel auxiliary subunit epsilon
Synonyms:	C1orf101; C10orf101
SwissProt:	Q5SY80
IHC positive control:	Human liver cancer

IHC Recommend dilution:

30-150

