

兔抗 CBLN1 多克隆抗体

中文名称: 兔抗 CBLN1 多克隆抗体

英文名称: Anti-CBLN1 rabbit polyclonal antibody

储 存: 冷冻 (-20℃)

抗 原: CBLN1

宿 主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a cerebellum-specific precursor protein, precerebellin, with similarity to the globular (non-collagen-like) domain of complement component C1qB. Precerebellin is processed to give rise to several derivatives, including the hexadecapeptide, cerebellin, which is highly enriched in postsynaptic structures of Purkinje cells. Cerebellin has also been found in human and rat adrenals, where it has been shown to enhance the secretory activity of this gland.
Applications:	ELISA, WB, IHC
Name of antibody:	CBLN1
Immunogen:	Synthetic peptide of human CBLN1
Full name:	cerebellin 1 precursor
SwissProt:	P23435

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	21 kDa
WB Positive control:	HUVEC cell lysate
WB Recommended dilution:	200-1000



