

兔抗 CDON 多克隆抗体

中文名称：兔抗 CDON 多克隆抗体

英文名称：Anti-CDON rabbit polyclonal antibody

别名：cell adhesion associated, oncogene regulated; CDO; CDON1; HPE11; ORCAM

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CDON

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a cell surface receptor that is a member of the immunoglobulin superfamily. The encoded protein contains three fibronectin type III domains and five immunoglobulin-like C2-type domains. This protein is a member of a cell-surface receptor complex that mediates cell-cell interactions between muscle precursor cells and positively regulates myogenesis.
Applications:	ELISA, IHC
Name of antibody:	CDON
Immunogen:	Fusion protein of human CDON
Full name:	cell adhesion associated, oncogene regulated

Synonyms:	CDO; CDON1; HPE11; ORCAM
SwissProt:	Q4KMG0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	25-100



