

兔抗 CNP 多克隆抗体

中文名称：兔抗 CNP 多克隆抗体

英文名称：Anti-CNP rabbit polyclonal antibody

别名：CNP1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CNP

反应种属：Human, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May participate in RNA metabolism in the myelinating cell, CNP is the third most abundant protein in central nervous system myelin.
Applications:	ELISA, WB, IHC
Name of antibody:	CNP
Immunogen:	Synthetic peptide of human CNP
Full name:	2',3'-cyclic nucleotide 3' phosphodiesterase
Synonyms:	CNP1
SwissProt:	P09543
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human gastric cancer
IHC Recommend dilution:	25-50

WB Predicted band size:	47 kDa
WB Positive control:	Human cerebrum tissue, Human cerebella tissue lysates
WB Recommended dilution:	500-2000



